

## CASKET TRAY with BAG SYSTEM



### Step # 1

With the casket tray on a flat surface, unfold the tray liner and lay it in the bottom of the tray.



### Step # 2

Place the casket in the tray and drape the casket bag over it and the tray.



### Step # 3

Leaving the excess bag hanging over the tray, gently tuck a small amount of the bag between the tray and the casket.

## CASKET TRAY with BAG SYSTEM page 2



### Step # 4

With the hooked side of the seal facing out, start in the middle of the tray and work the seal over the edge of the tray capturing the bag between the tray and the seal.



### Step # 5

Work your way around the tray, installing the seal until you get back to the starting point. You may need to slide the end forward or back to get the ends to meet.



### Step # 6

Once the seal is completely installed, carefully tuck the excess bag between the casket and the tray or trim the excess off leaving at least 2" below the seal.

Now place the crypt beads onto crypt floor, and slide the casket in.

### Casket Tray with Bag System \$125.00 each

Each Casket Bag System contains:

- 1 Polyethylene casket tray
- 1 absorbent liner
- 1 heavy 4mm bag with the check valve installed
- 1 rubber seal
- 1 one pound bag crypt beads

\*The Bag is available in either clear plastic (shown above) or Opaque White.

***If you are not completely satisfied with the Casket Bag System just return it for a refund!***

**C.A. Joseph Co. Inc.**  
13712 Old Fredericktown Rd.  
East Liverpool, Ohio 43920  
Phone: (330) 385-6869  
Fax: (330) 385-5038  
[www.cajoseph.com](http://www.cajoseph.com)